VIDA :

CRF15: Stool Accession (STA)

Web Version: 1.0; 1.04; 16Nov15

Specimen ID (SPEC_ID):

Is this specimen for VIDA or VIDA-Plus? (Select one)(STAVIDAS)

VIDA VIDA-Plus

Rectal swab sample collected?(STARECSW)

Whole stool sample collected?(STASTOOL)

No

Yes

Specimen receipt:

Time processed in laboratory (24 hour clock): (STASPTM) (hh:mm)

Date: (STASPDT) (ddMMMyyyy)

Properly labeled?(STALABEL)

Container tightly shut?(STACONST)

No
Yes

Temperature okay (cold pack frozen)?(STATEMP)

No
Yes

Sufficient weight/volume?(STAWVSF)

No Yes Actual weight: (STAWGT)

(xxx.xx) g

(minimum 4 g or ml for VIDA; 2 g or ml for VIDA-Plus)

Acceptable for accession? (STAACCES)

No Yes

Aliquoted for:

Viruses (ELISA) (STAVIRUS) No Parasites (STAPARA) No Yes Virus PCR (STAVRPCR) No

Yes

Yes

Archive 2(STAARCH2) No Archive 3(STAARCH3) No

Yes Yes

No

Specimen plated on:

Archive 1(STAARCH1)

SS/XLD (STASSXLD) No Yes TTGA/TCBS (STATTGAT) No Yes CAMPY (STACAMPY) No Yes

MaC(STAMAC) No Yes BP/APW(STABPAPW) No Yes ABA/CIN/RYAN(STAACR) No Yes

Technician: (STATECHN) (STATCNDT)

(ddMMMyyyy)

QC/Supervisor: (STASUPR)

Reviewed: (STASUPDT) (ddMMMyyyy)