VIDA

CRF17: Polymerase Reaction (PCR)

Web Version: 1.0; 1.04; 16Dec16

Specimen ID (SPEC_ID):

Date of Test: (PCRDT) (ddMMMyyyy)

Is this specimen for VIDA or VIDA-Plus? (Select one) (PCRVIDAS) VIDA VIDA-Plus

Multiplex PCR #1

(PCRM1TST) Test completed Not performed

estA (ST)(approx. 147 bp): (PCRM1ST)

eltB (LT)(approx. 508 bp): (PCRM1LT)

bfpA(approx. 300 bp): (PCRM1BFP)

aatA(approx. 630 bp): (PCRM1AAT)

positive

Negative

Negative

Negative

Negative

Negative

eae(approx. 881 bp): (PCRM1EAE) Positive Negative

eae(approx. oor bp). (PCRIVITEAE)

<u>Duplex PCR</u>

(PCRDPTST) Test completed Not performed

eae (approx. 384 bp): (PCRDPEAE)

ST (approx. 166 bp): (PCRDPST)

Positive Negative Negative

bfpA monoplex:

(PCRBFTST) Test completed Not performed

bfpA (approx. 300 bp): (PCRBFBFP) Positive Negative

Multiplex PCR #2

(PCRM2TST) Test completed Not performed

Stx-2 (approx. 443 bp): (PCRM2ST2)

eae (approx. 377 bp): (PCRM2EAE)

Sen (approx. 310 bp): (PCRM2SEN)

Stx-1 (approx. 220 bp): (PCRM2ST1)

Efa-1 (approx. 165 bp): (PCRM2EFA)

Positive Negative Negative

Technician: (PCRTECHN) (PCRTCNDT)

(ddMMMyyyy)

QC/Supervisor: (PCRSUPRV)

Reviewed: (PCRSUPDT) (ddMMMyyyy)